



CENTER FOR DISRUPTIVE
MUSCULOSKELETAL INNOVATIONS

Transporter Table System

Anand Agarwal, MD
David Dick
Vijay K. Goel, PhD
Adam MacMillan (IAB Mentor)
University of Toledo

This project funded by Spider Spine
Toledo, Ohio


WWW.NSFCDMI.ORG

This project aims to build a system for moving a patient between a gurney and operating room table without lifting.

- Protecting patients from handling injuries
- Protecting staff from Musculoskeletal Disorders (MSDs)
- Based on US Patent #9,463,127 has been granted to Leon Hochman, MD, Nancy Valvona, and Stephen Nicolato, West Bloomfield, MI.

¹<https://www.bls.gov/news.release/osh2.nr0.htm>

Progress – Computer Modelling



A solid model of one implementation has been made.

- This design consists of units made to mount on existing gurneys and beds
- Roller wheels extend and retract using air pressure
- Provision is made to lock the units together during the transfer

Progress – Prototype

The prototype has the following characteristics:

- Capable being mounted on an existing OR table and gurney
- Retractable wheel sets on the gurney and table
- A pallet to carry the patient between the tables
- Locks to prevent the pallet from moving when not in use
- A system to lock the gurney and table together during the transfer
- Means to direct the movement of the pallet to prevent skewing
- Continued access to existing gurney and table features such as accessory and guard rails.

Progress – Prototype

CENTER FOR DISRUPTIVE
MUSCULOSKELETAL
INNOVATIONS



The prototype system mounted on an OR table and gurney in the ECore Bio-Skills Lab. Here the “patient” is on the gurney and the rollers are retracted on both units


Progress – Prototype

CENTER FOR DISRUPTIVE
MUSCULOSKELETAL
INNOVATIONS



- Units latched together
- Air source connected
- Rollers extended

Progress – Prototype



Here the latch has been depressed, allowing the patient pallet to start moving towards the OR table

- The guide rails maintain the pallet alignment between the units.
- The latches engage the guide rails:
 - Prevent the pallet from overtravel
 - Secure the pallet during transport or surgery

Progress – Prototype

CENTER FOR DISRUPTIVE
MUSCULOSKELETAL
INNOVATIONS



Transfer complete

Progress – Prototype

CENTER FOR DISRUPTIVE
MUSCULOSKELETAL
INNOVATIONS



Evaluate the prototype in terms of ease of use, safety, and aesthetics:

- Demonstrate patient transfers at a design weight
- Measure the effort required to effect transfers
- Seek feedback from experienced staff